



# When to change?

Use until CO<sub>2</sub> exceeds  
5mmHg **or** FiCO<sub>2</sub> 0.5% vol.

**AMSORB® Plus** is **NOT**  
the same as sodalime.

No production of CO or Compound A (even when desiccated). Low and minimal flow anesthesia is considered **safer** with AMSORB® Plus\*.

Discard as non-hazardous waste.  
**Amsorb® Plus** exhausted pH <12.5, is **safe** for landfill and breaks down to harmless organic compounds.



\* Published papers available on request.



# When to change?

1. Remove packaging and install canister on anaesthesia machine.
2. Use until  $\text{FiCO}_2$  exceeds 0.5% or 5 mmHg  $\text{CO}_2$  breakthrough. **Colour will extend much further** through canister than with soda lime (see image above).
3. Discard as non-hazardous waste: **AMSORB® Plus** expired  $\text{pH} < 12.5$ , is safe for landfill and breaks down to harmless organic compounds.
4. MH (Malignant Hyperthermia) cases and Biohazard surgeries: **Follow existing hospital protocol**, including changing of the canister.

**Always** use capnometry to determine exhaustive state.



**Colourisation is normal** and expected and may reach 50% of the total absorber height prior to  $\text{CO}_2$  breakthrough.

[armstrongmedical.net](http://armstrongmedical.net)

