- INSTRUCTIONS FOR USE: 1. Connect each leg of WYE tubing to gas sources, providing flow ranges. 2. Connect Neo-Tee to WYE. 3. Adjust each flow device for desired flow range. NOTE: Flows from each device provide combined flow on outlet.

- WARNINGS: 1. Output of WYE flows are combined. 2. If using blended gases (i.e. air/O2), see chart below.

Table: Oxygen Flow Rate in LPM. Grid with flow rates for 0-15 LPM across 0-15 LPM.

Table: Iltflowhastighed i liter pr. minut. Grid with flow rates for 0-15 LPM across 0-15 LPM.

Table: Zuurstof stroomsnelheid in LPM (liters per minuut). Grid with flow rates for 0-15 LPM across 0-15 LPM.

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Table: Sauerstofffluss LPM (Liter pro Minute). Grid with flow rates for 0-15 LPM across 0-15 LPM.

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Table: Valore flusso ossigeno espresso in LPM (litri per minuto). Grid with flow rates for 0-15 LPM across 0-15 LPM.

NOTE: DATA IS LISTED AS "FLOW RATE IN LITERS PER MINUTE" L, "OXYGEN CONCENTRATION" % NOTE 2: OXYGEN CONCENTRATION IS BASED ON AIR CONTAINING 20.95% OXYGEN

- 3. Do not exceed 15 LPM total flow. 4. Always monitor output flows with oxygen analyzer, to confirm oxygen concentration. 5. To prevent flow restriction or back pressure, do not use additional connectors or swivels. 6. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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