SAFIRA

Instructions for Use (GB)
Instrucciones de Uso (ES)
Mode d'emploi (FR)

MEDOVATE
Developing Innovation

(English)**Definition:**

SAFIRA consists of a Sterile Syringe (20ml), Driver and Foot Pedal.
 The SAFIRA Syringe may be manually engaged and disengaged to the gearing of the Driver unit.
 The battery operated Driver is activated by means of a cable connected Foot Pedal.
 Flowrates are limited by design to a maximum of 0.5ml/sec.
 Approved needle types and extension tubing may be attached to the connector at the end of the SAFIRA syringe.
 SAFIRA is to only be used in conjunction with the single use SAFIRA Syringe .
 SAFIRAs internal battery power will perform 200 procedures.

Indications for Use:

The Medovate SAFIRA system is intended for use by trained clinicians to administer local anaesthetics below a specified pressure threshold to a target nerve bundle for regional anaesthesia.

Approved Needle Types:

Anaesthesia needles, which fall into the following ranges, have been approved for use with the SAFIRA system:

- The minimum sized needle gauge approved for use with this system is 22G
- The maximum length of the needle approved for use with this system is 120mm

Contraindications:

SAFIRA is not intended for the following uses:

- Intravascular delivery
- Delivery of blood, blood products, lipids, fat emulsions or Total Parenteral Nutrition (TPN)
- Infusion of fluids that will enter or contact circulatory blood or cerebrospinal fluid
- Delivery of life-supporting medications where under- or over-delivery may cause serious injury or death
- Use with Neonates and Paediatrics

Warnings:

1. Sterile technique should be used at all times during syringe filling, needle introduction and connection. Aseptic technique should be used for removal.
2. Medications or fluids must be administered per instructions provided by the drug manufacturer. Physician is responsible for prescribing drug based on each patient's clinical status (such as age, body weight, disease state of patient, concomitant medications, etc)

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MEDOVATE
 Developing Innovation

Medovate Ltd.
 The Workplace
 Camboro Business Park
 Girton
 Cambridgeshire
 United Kingdom
 CB3 0QH

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Warnings: (continued)

3. Make sure the medication being infused is approved for Regional Anaesthesia / PNB's (e.g. Lidocaine); follow all labelling instructions for medication use.
4. SAFIRA is MRI unsafe
5. Opening injection pressure varies among different tissues, being the highest when the needle tip is lodged in low compliance tissue (e.g., roots of brachial plexus, tendon) and lowest when injected in the soft connective tissues (e.g., adipose tissues, perineural space). SAFIRA limits the injection pressure to less than 20psi.
6. The SAFIRA Syringe is designed for single use in one patient and must not be re-sterilised, reconditioned or re-used. Re-use of the SAFIRA Syringe risks infection (due to accumulation of pathogens in the device that are subsequently injected) or potentially hazardous drug effects (due to residual drug compounds contaminating the intended drug and then being injected).
7. No modification of this equipment is allowed

Precautions:

FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN


1. SAFIRA has no user serviceable parts. **Do not attempt to repair or alter the device.**
2. The SAFIRA Syringe component is disposable and must be discarded after use in accordance with hospital, administrative and /or local government policy.
3. Do not use the SAFIRA Syringe if the packaging is open or damaged or if the sterile barrier is compromised.
4. The SAFIRA Syringe component is single use only.
5. The Driver component is a limited re-use device.
6. Flow rates may vary due to fill volume, viscosity and/or drug concentration, positioning the driver above or below the injection site, and temperature.
7. Start delivery within 8 hours of filling the SAFIRA Syringe . Storage of a filled SAFIRA Syringe component for more than 8 hours may result in slower flow rates.
8. SAFIRA is not made with natural rubber latex
9. SAFIRA is only to be operated by a trained Health Care Professional.
10. The SAFIRA Syringe is not intended for measurement use.

SAFIRA Instructions For Use:


- 1. SAFIRA consists of 3 separate components. They are the **SAFIRA Syringe**; the **Driver**; and, the **Foot Pedal** (all 3 are pictured and labelled below).

The SAFIRA Syringe (Luer connection variant shown, NRFit connection variant is also available)

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



SAFIRA Syringe Definition

The SAFIRA Syringe is available in two variants defined by it's barrel connection interface. One variant uses a Luer connection interface and the other an NRFit connection interface. All variants of the SAFIRA Syringe are compatible with the SAFIRA Driver. The different variants can be identified in the following way:

SAFIRA Syringe (Luer Connection)

SAFIRA Syringe (NRFit Connection)

- 2. Open the sterile package and remove SAFIRA Syringe. Use normal hospital technique to fill and prepare the syringe.
- 3. Connect the appropriate needle to the fitting of the syringe. Ensure that the needle is connected securely to the syringe. Once secured, prime the needle and tubing set to remove the air.


Diagram #1: Needle prepared to connect to fitting of SAFIRA syringe

- 4. To attach the SAFIRA Syringe to the Driver; **remove protective foam cap from the Driver gear prior to connecting (see below)**. Align SAFIRA Syringe with Driver as shown in Diagram #2. Once aligned, gently press the SAFIRA Syringe into the Driver. It will be seated correctly when you hear it “click” into place.


Diagram #2: SAFIRA ready for the SAFIRA Syringe

- 5. Attach the Foot Pedal to the Driver. The plug end of the cord goes into the round receptacle opening as shown in Diagram #3. Place Foot Pedal on the floor.


Diagram #3: Foot Pedal plugged in

- 6. Press the infuse part (green end, furthest from the cable) of the Foot Pedal to check that the tubing is correctly primed and the SAFIRA system is operating as expected.

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7. Place the Foot Pedal on the floor below the surgical site in a convenient location for the Anaesthesiologist. See Diagram #4 for reference.


Diagram #4: Foot Pedal placed on floor below surgical site

8. Removal: In order to remove the SAFIRA Syringe from the Driver, simply press the square teal button on the side of the Driver as shown in Diagram #5. The SAFIRA Syringe will immediately release from the Driver.


Diagram #5: Press square release button to remove Syringe

9. The Foot Pedal has two colours: Green = Infusion; and, Yellow = Aspiration as shown in Diagram #6.


Diagram #6: Foot Pedal showing Infusion and Aspiration functions

10. **Driver Indicator Lights:** On the top edge of the Driver (refer to **Diagram #7**) colour indicator lights are displayed representing various functions in action or functions that require additional action. The indicators are as follows:

LED Colour	State	Meaning	Operator Response
Red	Flashing (Infusion has stopped)	Warning: <u>Immediate response</u> by operator is required	Operator: Infusion has stopped; 1. Verify needle patency and/or reposition needle to allow low pressure infusion 2. Reset SAFIRA (see steps 12 and 13 for details)
Red	Steady (End of Life)	Warning: <u>Immediate response</u> by operator is required.	Operator: Immediately disconnect the syringe and either replace Driver or continue the procedure manually.
Green	Steady (Infusion Function Active)	Normal: <u>System operating as designed indicating infusion activity</u>	Operator: Continue infusion activity until medical result is achieved.
Yellow	Steady (Aspiration Function Active)	Caution: <u>Aspiration occurring</u> - system operating as intended	Operator: Follow Standard Precautions during use of aspiration function
No Active Light	N/A	Normal: <u>System is ready for infusion or aspiration</u>	Operator: N/A


Diagram #7: Top edge of driver, illustrated with all indicators.

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11. Begin the procedure following standard hospital guidelines. Throughout the injection SAFIRA is designed to limit the infusion pressure up to a maximum of 20 psi.
12. Should an infusion stoppage be encountered during the procedure, it may indicate an intraneural injection or a blockage in the fluid path. Verify needle patency and/or reposition needle to allow low pressure infusion and then **reset the Flashing Red Light/Infusion Stopped** on the Driver.
13. To **reset the Flashing Red Light / Infusion Stopped** on the Driver, follow either of the two (2) steps indicated below:
 - a. **Aspirate** (Yellow side closest to the cable of the Foot Pedal) until red lights go out: or
 - b. **Release Sterile Syringe** (see **Diagram #5**) and re-click into place

Now the procedure can be continued.
14. When the procedure is complete, follow standard hospital protocol for the following:
 - a. Dispose of the Sterile Syringe using acceptable standard practice for biohazard waste
 - b. Driver and Foot Pedal components should be stored according to hospital practice for reuse
15. **End of Life:** Once the end of life indicator on the Driver (steady red lights) comes on, the component should be disposed of in accordance with hospital procedure. Do not start a procedure if the end of life indicator is already on from the beginning. If the end of life indicator lights up during a procedure, immediately disconnect the syringe and either replace the Driver or complete the procedure manually.

Disposal:

For SAFIRA components, follow standard hospital protocol for disposal.

1. **SAFIRA Syringe:** Use acceptable practice for biohazard waste.
2. **Driver:** Should not be disposed of as unsorted municipal waste. Dispose of unit in line with local guidelines.
3. **Foot Pedal:** Should not be disposed of as unsorted municipal waste. Dispose of unit in line with local guidelines.

Reuse:

The SAFIRA **Syringe** component **is not reusable** and must be discarded using standard biohazard disposal procedures. The SAFIRA **Driver** component **is reusable** and should be wiped with gauze soaked in at least 70% isopropyl alcohol. The SAFIRA **Foot Pedal** component **is reusable** and should be wiped with gauze soaked in at least 70% isopropyl alcohol. Failure to follow hospital guidelines may result in patient cross contamination.

How Supplied:

The SAFIRA Syringe component is provided sterile; the contents are sterile unless the package is opened, damaged or the expiry date has been exceeded.

Use Environment:

SAFIRA is designed for use by a physician, in either a hospital or surgical centre environment. The device is not intended to be used outside of stated environments.

Environmental Operating Conditions, Transport and Storage Between Uses:

(Driver and Foot Pedal)

Temperature Range (Product Use): 50F (10C) to 104F (40C)

Temperature Range (Shipping & Handling): 14F (-10C) to 104F (40C)

Humidity Range: 10% to 95%, noncondensing

Atmospheric Pressure: 500 to 1060 Millibars

Device Type:

SAFIRA is a Type BF device.

The SAFIRA components are not conductive and can be immediately released from patient. The needle and tubing (which is not supplied by Medovate) attached to the SAFIRA is the part in physical contact with the patient and can be immediately released from the patient.

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Transient.

Power Supply:

Two (2) AAA 1.5V alkaline batteries power the Driver component.

Trouble Shooting:

1. See **Diagram #7** on Page 7 describing actions required should either flashing red light or steady red light appear.
2. **If unit stops working during the middle of a procedure, simply press the teal button on the side of the Driver component and release the SAFIRA Syringe. You may complete the procedure using the SAFIRA Syringe in the traditional manual mode.**
3. Should the Driver component fail to operate either before or during a procedure, please return to Medovate for replacement. **The device should not be taken apart or repaired by anyone other than authorized Medovate personnel.**
4. Should a serious incident occur during the use of or as a result of the SAFIRA system, please contact Medovate.
5. See Medovate website for details of local Medovate representative.

End of Section

Important Technical Information

Electromagnetic Compatibility (EMC):

This product needs precautions regarding EMC and needs to be installed into service according to the EMC information provided. This unit can be affected by portable and mobile RF communication equipment.

1. Do not use a mobile phone or other devices that emit electromagnetic fields near the unit. They may result in incorrect operation of the unit.
2. **CAUTION:** SAFIRA should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this unit should be observed to verify normal operation in the configuration in which it will be used.

Electro Magnetic Emission:

Guidance and Manufacture's Declaration: Electromagnetic Emission		
SAFIRA is intended for use in the electromagnetic environment specified below. The customer and/or user of SAFIRA should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic Environment: Guidance
RF Emissions CISPR 11	Group 1	SAFIRA uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Group A	SAFIRA is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and Manufacture's Declaration: Electromagnetic Immunity			
SAFIRA is intended for use in the electromagnetic environment specified below. The customer and/or user of SAFIRA should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment: Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±2, 4, 8, 15 kV air	±8 kV contact ±2, 4, 8, 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power Frequency (50Hz/60Hz) Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U _T is the a.c. mains voltage prior to application of the test level.			

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Electro Magnetic Immunity: (continued)

Guidance and Manufacture's Declaration: Electromagnetic Immunity

SAFIRA is intended for use in the electromagnetic environment specified below. The customer and/or user of SAFIRA should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment: Guidance
Conducted RF IEC 61000 -4-6	3V rms (1kHz - 80%) 150kHz – 80MHz 6V rms ISM and amateur radio bands	3Vms 6V rms	Portable and mobile RF communication equipment should be no closer to any part of SAFIRA, including cables, than the recommended separation distance of 300mm. Field strength from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b .
Radiated RF IEC 61000-4-3	3V/m (1kHz - 80%) 80MHz – 2.7GHz	3V/m	Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: The guidelines may not apply in all situations, Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Electro Magnetic Immunity: (continued)

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which SAFIRA is used exceeds the applicable RF compliance level above, SAFIRA should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating SAFIRA.

^b Over the frequency range of 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

End of Section

Guide to SAFIRA References:


- 900016 (Component Level—stands for the SAFIRA Luer Syringe component)
- 900058 (Component level—stands for the SAFIRA NRFit Syringe component)
- 900029 (Component Level—stands for the Driver component)
- 900044 (Component Level—stands for the Foot Pedal Component)

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


Pédale de pied


Guide to Symbols:

#	Symbol	Symbol Meaning
1		Component reference number
2		Lot number or batch code
3		Do not use if package is damaged
4		Single use only
5		Use by date
6		Date of manufacture
7		Caution: US Federal Law restricts this device to sale by or on the order of a physician
8		Consult Instructions for Use (IFU) (Recommended)
9		Product sterilised using ethylene oxide
10		Consult Operating Instructions or IFU (Mandatory)
11		Type BF: Provides protection against electrical shock
12		Electromagnetic emissions

Guide to Symbols Continued:

#	Symbol	Symbol Meaning
13		Temperature limitation
14		Do not re-sterilise
15		MRI Unsafe
16		Manufacturer
17		WEEE compliant
18		RoHS compliant
19		CE Mark
20		Luer Syringe
21		NRFit Syringe

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Warranty Statement:

The product and each component of its system (hereinafter "the product") have been designed, manufactured, tested and packaged with all reasonable care. However, Medovate has no control over the conditions under which the product is used and a disturbance of the intended function of the product may occur for various reasons. In this respect, the warnings in the product publication/instructions for use are expressly to be considered as an integral part of this disclaimer and provide more detailed information. For this reason, Medovate disclaims all warranties, expressed or implied regarding the product, including but not limited to, any warranty of merchantability or fitness for a particular purpose of the product. Product descriptions or user guidelines in publications do not constitute any expressed representation, or any expressed or implied warranty. Medovate is not liable for any direct, incidental or consequential damages or medical expenses caused by any use, defect, failure or malfunction of the product whether the claim is based on contract, warranty, tort or otherwise. This does not apply in the case of intention or in the case of gross negligence of legal representatives or executive staff of Medovate. In commercial transactions relating to merchants, the liability is limited to the compensation of typical damages; compensation for any untypical or incidental damage is excluded. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable in the respective country. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Medovate's legitimate interest in limiting liability or warranty without infringing any mandatory provisions of applicable law. No person has any authority to bind Medovate to any warranty or liability regarding this product.

SAFIRA...SAFer Injection for Regional Anaesthesia

The SAFIRA is manufactured for:

MEDOVATE
Developing Innovation

Medovate Ltd.
The Workplace
Camboro Business Park
Girton
Cambridgeshire
United Kingdom
CB3 0QH

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