


ONLY ONE



NIFometer
(to -40 cm H₂O)
10-340-00
20/Box



NIFometer
(to -60 cm H₂O)
10-340-02
20/Box

Single-Patient-Use

NIFometer®

NEGATIVE INSPIRATORY FORCE METER

The quick & convenient way to measure NIF.

QUICK AND EASY

- Small, compact, self-contained. Connects directly to ET tube eliminating need for separate tubing and connectors.
- Includes built-in occlusion button that you simply depress during test.

SINGLE PATIENT USE

- Helps eliminate risk of cross-contamination.

ACCURATE

- Accurate to ± 2 cm H₂O -20 to -30 cm H₂O (NIF 0-40) & ± 3 up to -15 cm H₂O & ± 5 from -15 to -60 cm H₂O (NIF 0-60).
- No calibration required.
- Includes memory indicator pointer that records highest force reached and can then be reset.

CONVENIENT

- Lightweight weighing only 30 Grams.
- Provides unobstructed view of patient during test.
- Small and compact, fits in pocket – readily available for added efficiency.
- Also available in kit with elbow, mouthpiece and nose clip.

ONLY ONE

Single-Patient-Use

NIFometer®

NEGATIVE INSPIRATORY FORCE METER

The quick & convenient way to measure NIF.



NIFometer
(to -40 cm H₂O)

NIFometer
(to -60 cm H₂O)



NIFometer Kit

Available with NIFometer
to -40 cm H₂O & to -60 cm H₂O

SPECIFICATIONS

| | |
|----------|---|
| Range | To -40 cm H ₂ O |
| Variance | ± 2 cm H ₂ O -20 to -30 cm H ₂ O |
| Fitting | 15 mm I.D./22 mm O.D. |
| Weight | 30 grams nominal |
| Range | To -60 cm H ₂ O |
| Variance | ± 3 up to -15 cm H ₂ O & ± 5 from -15 to -60 cm H ₂ O |
| Fitting | 15 mm I.D./22 mm O.D. |
| Weight | 30 grams nominal |

ORDERING INFORMATION

| Part # | Description | Packaged |
|-----------|---|----------|
| 10-340-00 | NIFometer (to -40 cm H ₂ O) Single-Patient-Use | 20/Box |
| 10-340-01 | NIFometer Kit, (to -40 cm H ₂ O) Includes: elbow, mouthpiece & nose clip | 10/Box |
| 10-340-02 | NIFometer (to -60 cm H ₂ O) Single-Patient-Use | 20/Box |
| 10-340-03 | NIFometer Kit, (to -60 cm H ₂ O) Includes: elbow, mouthpiece & nose clip | 10/Box |

Patent # US 6,854,334 B2