

NEVER STOP A PROCEDURE AGAIN TO CHANGE AN EMPTY CO₂ TANK


Switch from your insufflator's empty cylinder to a full cylinder in *seconds*.

ADVANTAGES OF THE INSUFFLATOR DUAL TANK KIT

- Never stop a procedure to change CO₂ tanks again!
- Connect two CO₂ tanks simultaneously to any insufflator.
- Easy-to-operate handle to switch CO₂ sources.
- Low profile, lightweight set-up.
- Dual tank kits are configured to fit any brand of insufflator applications.


Mercury has kits available to accommodate all insufflator and tower combinations.

Contact Mercury Medical's Customer Service at
(800) 835-6633 for additional information.



BEFORE

AFTER



*“switch-over”
in seconds*